


Chlamydia CEL IFA

The Cellabs Chlamydia CEL IFA kit is a direct immunofluorescence assay for the detection of *Chlamydia trachomatis* organisms in patient specimens or in culture. It uses a monoclonal antibody specific for an outer membrane protein (OMP) of *Chlamydia trachomatis*.



Superior clinical reliability

100% NPV and 100% PPV vs culture method (n=270)
100% NPV and 94% PPV vs culture & Giemsa (n=400)

Most reliable and affordable

- No culture required, direct observation
- Cost effective alternative to PCR methods
- Simple and rapid (30 minutes incubation)



Cellabs Pty Ltd

P.O. Box 421 Brookvale
NSW 2100 Australia
Ph: + 61 2 9905 0133
Fax: + 61 2 9905 6426
Website: www.cellabs.com.au
Email: sales@cellabs.com.au