


# Giardia CEL IFA

Giardia-CEL is a direct IFA kit using a monoclonal antibody specific for a cell wall protein of *Giardia* cysts in faecal specimens, water or sewage samples.


## Clinically reliable

Reliable *Giardia* confirmatory test  
Control slide included in the kit  
Convenient for use on single or few specimens (no reagent waste)



## Ideal for environmental tests

Can be used on filtered water samples or sewage samples  
Excellent fluorescence and clearly observable morphology for easy pathogen detection



## Reliable *Giardia* confirmation

- Highly specific (96%) for *Giardia*
- Rapid results (30 minutes incubation)
- Ready-to-use reagents
- 50 tests kit format with long shelf life



### Cellabs Pty Ltd

P.O. Box 421 Brookvale  
NSW 2100 Australia  
Ph: + 61 2 9905 0133  
Fax: + 61 2 9905 6426  
Website: [www.cellabs.com.au](http://www.cellabs.com.au)  
Email: [sales@cellabs.com.au](mailto:sales@cellabs.com.au)